

Insemination Wagon DATASHEET JUNE 2017

Insemination wagon makes efficient artificial insemination works.
The wagon includes: 1 incubator, 1 tool box, 1 tissue box and 1 plastic drum.
There are 4-wheel and 2-wheel types for selection.



Artificial Insemination Wagon

- ①The incubator is used for storage of seminal fluid of pigs.
- ②The white drum is used for storage of seminiferous duct.
- ③White tissue box with zigzag open mouth for breaking the tissue easily.
- ④The storage box is for other tools.

Size: 1220*680*920mm
Incubator capacity:26L